

Brewer

Ergonomic Medical Seating

Designed to be a comfortable
fit for your team and
your practice.



Millennium Series

Ergonomically designed to meet your practice's needs.

Seating tends to be one of the most under appreciated aspects of a clinic. Yet it can dramatically impact the productivity and health of your entire staff. That's why our ergonomic experts work closely with clinicians to develop specific comfort solutions for their specialized requirements. The Millennium Series offers the thoughtful design and premium features that help make every day on the job seem less like work. Plus, a standard five-year warranty means added peace of mind.

Millennium Series

Premium specialty seating features seamless upholstery, an 18.5" (47 cm) wide, contoured seat cushion and adjustable contoured backrest unless otherwise noted. For stability, all incorporate a 23" (58.4 cm) five-leg, powder-coated black aluminum base and 2" (5.1 cm) dual-wheel nylon casters with polyurethane tread unless otherwise noted. Weight capacity: 250 lbs. (113 kg).

Task Models

Fully adjustable forward seat tilt and synchronized seat and backrest tilt support neutral posture to reduce fatigue and back pain. Pneumatic height adjustment: 18"-24.75" (45.7 cm-62.9 cm).

21435BA - w/ arm rests

21435B - w/o arm rests



Model 21435BA in Azure Blue

A - Arm rests B - Backrest



Lab Model 21521B in Adriatic upholstery

Surgeon Model 21340 in Black upholstery

Task Model 21435BA in Azure Blue upholstery

Lab Models

Features 18" (48.7 cm) foot ring standard and pneumatic height adjustment:
 – w/ glides: 20.75"-31" (52.7 cm-78.7 cm)
 – w/ casters: 21.75"-32" (55.2 cm-81.3 cm)

21520B - w/ glides
 21521B - w/ casters



Model 21521B in Adriatic

Surgical Models

Offers hands-free, foot-operated pneumatic adjustment:
 19"-24" (48.3 cm-61 cm)

21340B - 18.5" (47 cm) seat
 21340BV - 18.5" (47 cm) vacuum-formed seat
 21340 - 16" (40.6 cm) seat; seamless upholstery not standard
 21340V - 16" (40.6 cm) vacuum formed seat



Model 21340 in Black



Model 21340B in Wedgewood

Accessories

2102265 - set of 5 pressurized locking casters (sit to lock casters)
 3102277 - replacement casters (set of 5)

B - Backrest V - Vacuum formed/Seamless

Century Series

An ideal blend of quality, comfort and value.

Brewer's ability to work quality craftsmanship and comfort-driven features into affordable seating solutions has made our Century Series and Value Plus Series popular choices at medical facilities around the world.

Century Series

Offers the comfort of a 16" (40.6 cm) diameter, 4" (10.2 cm) thick poly-foam seat cushion with the dependability of a five-leg, 22.25" (56.5 cm) brushed aluminum accent

base. 20" - 27.5" (50.8 cm - 69.9 cm)

Adjustable Height

Models

- 11001BFR
- 11001BVFR
- 11001BFR-C133
- 11001FR
- 11001VFR
- 11001FR-C133



Model 11001BFR in Cabernet

Model 11001FR in Cabernet

17.75" - 24.5" (45.1 cm - 62.2 cm)

Adjustable Height

Models

- 11001B
- 11001BV
- 11001BLCD
- 11001BG*
- 11001BFR-C133
- 11001B-UL
- 11001
- 11001V
- 11001LCD
- 11001G*
- 11001-C133
- 11001-UL



Model 11001B in Black

Model 11001 in Black

* Height range: 16.75" - 23.5" (42.5 cm - 59.7 cm)

B - Backrest FR - Footring V - Vacuum formed/Seamless
 LCD - Locking casters G - Glides C133 - CAL-133 upholstery
 UL - Ultraleather™ upholstery



Model 11001BFR in Cabernet upholstery

Model 11001 in Black Satin upholstery

Value Plus Series

Designed for affordability, each features a 15" (38.1 cm) wide, 4" (10.2 cm) thick poly-foam seat cushion and five-leg, 22" (55.9 cm) black ABS plastic composite base.

17" - 21.25" (43.2 cm - 54 cm)
 Pneumatic Height Adjustment

Models

- 22500B
- 22500BV
- 22500B-C133
- 22500
- 22500V
- 22500-C133



Model 22500B in Ash

Model 22500 in Ash

17" - 21.25" (43.2 cm - 54 cm)
 Spinlift Height Adjustment

Models

- 22400B
- 22400BV
- 22400
- 22400V



Model 22400B in Ash

Model 22400 in Ash

B - Backrest FR - Footring V - Vacuum formed/Seamless
 LCD - Locking casters G - Glides C133 - CAL-133 upholstery
 UL - Ultraleather™ upholstery

Additional Practice Seating

Ergonomic solutions for a wide range of needs.

Round out your clinic seating needs with rugged, reliable designs from Brewer. It's a choice your staff and your patients will appreciate.

Tradition Series

This economical exam stool is made to stand the test of time. The Tradition Series features spinlift height adjustment from 19.5"-27.5" (49.5 cm-69.9 cm), a sturdy four-leg chrome base, 14" (35.6 cm) wide, 3.5" (8.9 cm) thick poly-foam seat cushion and 2" (5.1 cm) rubber wheel, ball bearing casters.

Model 23051

Side Chair Series

Perfect for exam room, front office or waiting room environments, these durable chairs feature steel tube frames and kick-out rear leg design to help prevent wall dings and dents. Their contoured backrests and 2" (5.1 cm) thick poly-foam seat cushions offer plenty of comfort. Plus, they meet BIFMA (Business and Institutional Furniture manufacturer's Association) Standards.

Models

1200 - w/ sturdy plastic arm rests

1250 - w/o arm rests



Model 1250 Deep Sapphire upholstery



Model 23051 in Arctic



Model 1200 in Ivy

Blood Drawing Chair Series

Seat alignment positions patients securely on the chair, while innovative angled armrests and easy-to-use knob height adjustments provide optimal arm position for each patient's blood draw. Features include two standard add-on work surfaces (to accommodate prep materials, paperwork and vials). Chair also meets BIFMA Standards.

An optional drawer (Model #101059) can be attached on either side of chair to keep prep materials nearby yet out of sight.

Model 1500

Accessories

101059 - Drawer option



Model 1500 in Cabernet

Upholstery Color Options

Standard upholstery color selections


For Brewer Side Chairs, Blood Drawing Chairs and Exam Stools



California Technical Bulletin 133 Upholstery Color Selections (CAL-133)

For Brewer Century Series Stools and Value Plus Series Pneumatic Stools


CAL TB 133 is a written legal standard outlining flammability test procedures for upholstered furniture. These vinyl color selections contain a flame blocking barrier that helps to prevent ignition, reflects heat and is flame retardant.



Standard upholstery color selections

For Brewer Side Chairs, Blood Drawing Chairs and Exam Stools

All Premium Plush Upholstery Colors are Ultrafabrics ULTRALEATHER™.



Colors may vary slightly depending on printing setting. For actual material swatches call 1•888•Brewer•1

The Power to Advance



Brewer seating solutions are proudly made in the USA. Our extensive line of solutions is manufactured to our rigorous quality standards and backed by best-in-class standard warranties. It's all part of our commitment to provide you with The Power to Advance.

For more information about the many ways we deliver The Power to Advance, contact your Brewer representative today.

Brewer

The Brewer Company, LLC
N88 W13901 Main Street, Suite 100
Menomonee Falls, WI 53051

1.888.273.9371
www.brewercompany.com

