



# VISTOR PRO LED

EXAM AND PROCEDURE LIGHT






BRINGING  
**VISION**  
TO THE OPERATING ROOM



PRODUCT  
SPECIFICATIONS

# LED TECHNOLOGY IS NOW OUTSIDE OF THE OPERATING ROOM

## CONFIGURATION OPTIONS



Available in single or dual ceiling mount, wall mount or mobile stand

The VISTOR PRO LED is another addition to the growing line of NUVO LED Lights. The state-of-the-art design provides unparalleled optics for ambulatory surgery centers, procedure rooms, trauma rooms, clinics and physician offices.

Control on/off and intensity settings from sterile handle.



The VISTOR PRO LED provides superior illumination and unmatched shadow control using 36 high efficiency ultracool LED's.

## PRODUCT SPECIFICATIONS

<b>Light Intensity</b>	64,500 lux
<b>Color Temperature</b>	4500 K
<b>Intensity Settings</b>	5(50% – 100%)
<b>Controls in Sterile Field</b>	Yes
<b>LED Life Efficiency</b>	> 30,000 hrs.



*A member of the  
Medical Illumination International Group*

5368 Kuhl Road • Erie, PA 16510-4703

(800) 663-1152 • (814) 899-4220

Fax: (814) 899-1410

[www.nuvosurgical.com](http://www.nuvosurgical.com)

