



NUMERICAL INDEX

022.....	21	8642.....	8
055.....	13 & 21	8648.....	9
066.....	13 & 21	8672.....	10
076.....	21	8711.....	23
076FF.....	13	8711-A.....	23
089.....	28	8720.....	23
24P.....	12 & 20	8720-A.....	23
36P.....	12 & 20	8805.....	25
42P.....	12 & 20	8805-A.....	26
48P.....	12 & 20	8805-AF.....	26
72P.....	12 & 20	8821.....	25
36Q.....	12	8821-A.....	26
42Q.....	12	8821-AF.....	26
48Q.....	12	8822.....	25
72Q.....	12	8822-A.....	26
8024.....	16	8822-AF.....	26
8024-99.....	18	8900-A.....	28
8036.....	16	8920.....	28
8042.....	17	8921.....	24
8042-99.....	19	8921-A.....	27
8048.....	17	8921-AF.....	27
8048-99.....	19	8922.....	24
8072.....	18	8922-A.....	27
8224.....	16 & 18	8922-AF.....	27
8236.....	16	8923-A.....	28
8242.....	17	8950.....	24
8242-9.....	19	8950-A.....	27
8248.....	17 & 19	8950-AF.....	27
8272.....	18	9800.....	31
8324.....	6	9810.....	31
8336.....	6	9820.....	31
8342.....	8		
8348.....	9		
8372.....	10		
8624.....	6		
8636.....	6		

COLOR SELECTORS:  
In an effort to stay current, all color selections for individual products can be found online at: [www.clinton-ind.com](http://www.clinton-ind.com)



WHY CLINTON?

Buyers of medical equipment are faced with evaluating clinical needs and justifying expenditures on nearly every purchase. Now, Clinton makes it a whole lot easier on everybody with products that maximize performance and aesthetics while minimizing costs.

**Performance.** With nearly forty years in medical manufacturing, Clinton has met the evolving challenges of the industry with products that fit the needs of providers, while at the same time maintaining a track record of durability that boasts the best warranty in the business.

**Aesthetics.** Nobody likes outdated anything, so Clinton has changed its products over the years to reflect the advances in medical practice, design and contemporary decor, making products with appeal in both form and function.

**Cost.** Bells and whistles are great, but they also drive up costs. On the other hand, why pay more if you don't have to? Clinton products fill needs instead of empty pockets, providing medical professionals with just the right products—at just the right prices.



**IMPORTANT NOTE TO DISTRIBUTORS DOING BUSINESS IN CALIFORNIA:**  
Clinton Industries products drop shipping or to our knowledge being reshipped into California will carry appropriate Safe Harbor Warnings in compliance with California Proposition 65.  
For more information on California Proposition 65, go to: [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov)





FASHION FINISH CASEWORK



Fashion Finish Casework by Clinton features great looks, premium performance and a tough, heat-sealed thermofoil surface, making Fashion Finish easy to clean and even easier on budgets.

*Fashion Finish Casework features:*

- Textured thermofoil doors, drawer fronts and sides
- MDF core construction
- Easy-clean, durable, surfaces
- Soft-close, door hinges
- 75 lb. capacity drawers
- Easy installation
- Select postform or quartz countertops
- Sinks and faucets are standard
- Available in select sizes



Choose from five contemporary colors specifically selected to span across the spectrum of modern medical design and decor.

**Sunlight Oak.**

Breathe in the airy feel of this light woodgrain surface and brighten your space with a contemporary blend of warm tones.

**Metropolis Gray.**

Gray is trending for a number of reasons. It blends perfectly with modern medical decor and adds a soft elegance wherever it is used.

**Chestnut Hill**

A classic blend of rich earthtones highlights this new style of woodgrain. Chestnut Hill, like its namesake, adds a touch of class and comfort to your professional space.

**Arctic White**

Bright with crisp lines and a smooth, snow white surface, Arctic White is the ultimate in the “clean room” look, or use as a sharp contrast to a dark accent wall.

**Twilight.**

At sunset, light yields to the dark tones that are found in the modern feel of our Twilight woodgrain.

Use Twilight as a bold contrast to lighter palette decor.

**NOTE:** Fashion Finish Casework does not include a countertop. Countertops must be ordered separately.







Wall Cabinet with 1 Door

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Easily mounts with hook and wall rail mounting system
- Shipped assembled for easy installation

8324	width	depth	height
	24"	12"	24"
	61 cm	30 cm	61 cm

Base Cabinet with 2 Doors & 2 Drawers

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Select postform laminate or quartz countertops
- Shown with Postform top and top mount sink and faucet (8624 *not available with quartz top*)
- Shipped assembled for easy installation

8624	width	depth	height
	24"	18"	35"
	61 cm	46 cm	89 cm



Wall Cabinet with 2 Doors

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Easily mounts with hook and wall rail mounting system
- Shipped assembled for easy installation

8336	width	depth	height
	36"	12"	24"
	91 cm	30 cm	61 cm

Base Cabinet with 2 Doors & 2 Drawers

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Optional 5-drawers may be added to replace 2-drawers and door section
- Select postform laminate or quartz countertops
- Shown with optional quartz top and undermount sink and faucet
- Shipped assembled for easy installation

8636	width	depth	height
	36"	18"	35"
	91 cm	46 cm	89 cm







Wall Cabinet with 2 Doors

- 1 adjustable shelf
- Shipped assembled for easy installation
- Concealed, self-closing, adjustable, door hinges
- Easily mounts with hook rail mounting system

8342	width	depth	height
	42"	12"	24"
	107 cm	30 cm	61 cm

Base Cabinet with 3 Doors & 2 Drawers

- 2 adjustable shelves
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Optional 5-drawers may be added to replace 2-drawers and door section
- Select postform laminate or quartz countertops
- Shown with postform top and top-mount sink and faucet
- Shipped assembled for easy installation

8642	width	depth	height
	42"	18"	35"
	107 cm	46 cm	89 cm



Wall Cabinet with 2 Doors

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Easily mounts with hook and wall rail mounting system
- Shipped assembled for easy installation

8348	width	depth	height
	48"	12"	24"
	122 cm	30 cm	61 cm

Base Cabinet with 2 Doors & 2 Drawers

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Optional 5-drawers may be added to replace 2-drawers and 2-doors section
- Select postform laminate or quartz countertops
- Shown with postform top and top-mount sink and faucet
- Shipped assembled for easy installation

8648	width	depth	height
	48"	18"	35"
	122 cm	46 cm	89 cm



Wall Cabinet with 3 Doors

- 2 adjustable shelves
- Shipped in 2 units for easy installation
- All holes are pre-drilled and all hardware is provided to connect the 2 cabinets
- Concealed, adjustable, soft-close, door hinges
- Easily mounts with hook and rail mounting system

8372	width	depth	height
	72"	12"	24"
	183 cm	30 cm	61 cm

Base Cabinet with 6 Doors & 3 Drawers

- 3 adjustable shelves
- Shipped assembled for easy installation
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Optional drawers may be added or repositioned to customize your layout
- Select postform laminate or quartz countertops
- Shown with postform top and top-mount sink and faucet
- Shipped assembled for easy installation

8672	width	depth	height
	72"	18"	35"
	183 cm	46 cm	89 cm







**POSTFORM LAMINATE COUNTERTOP**

Postform laminate countertops offer a distinctive one-piece countertop that eliminates the seam between countertop surface, backsplash and front edge. Postform countertops look great, are economical, durable and easy to clean.

*Note: All postform tops come standard with top mount, stainless steel sink and dual lever chrome faucet. All tops without sinks are special order only. All 72" postform tops with sinks positioned other than centered are special order. 72" cabinet's tall doors MUST be positioned under sink to accommodate plumbing.*



**HANSTONE® QUARTZ COUNTERTOP**

HanStone is nonporous, has no edge seams, is twice as strong as granite and will not absorb liquids like other surfaces. It is heat, bacteria and stain resistant and is GREENGUARD Indoor Air Quality Certified®.

A 3 1/4" (8 cm) backsplash is included with countertop and is installed in the field.

*Note: All quartz tops come standard with undermount, stainless steel sink and single lever stainless steel faucet. All quartz tops without sinks are special order only. All 72" quartz tops with sinks positioned other than centered are special order. 72" cabinet's tall doors MUST be positioned under sink to accommodate plumbing. Quartz top not available for 24" cabinets.*

24P	24" Top with sink & faucet
36P	36" Top with sink & faucet
42P	42" Top with sink & faucet
48P	48" Top with sink & faucet
72P	72" Top with sink & faucet



WC	White Cararra
BA	Black Alicante
ME	Meteorite

36Q	36" Top with sink & faucet
42Q	42" Top with sink & faucet
48Q	48" Top with sink & faucet
72Q	72" Top with sink & faucet



IW	Ivory Wave
BC	Black Coral
RB	Royale Blanc

**HanStone™**  
Fine Quartz Surfaces

**TOP MOUNT STAINLESS STEEL SINK AND CHROME FAUCET (POSTFORM ONLY)**

A stylish, easy clean, top mounted, stainless steel sink and chrome faucet is standard with all postform countertops. Faucet is 14.5" high chrome with a 3" radius and dual, wing levers .

**Sink & Faucet included with Postform Top**

inside sink dimensions

10" x 12" x 6"

25 cm x 30 cm x 15 cm

18 ga stainless steel

**UNDERMOUNT STAINLESS STEEL SINK AND FAUCET (QUARTZ ONLY)**

Undermount sink easily installs under quartz top and sealed with appropriate sealer. Faucet is 14.5" high stainless steel, single lever with a 7.9" radius.


**Sink & Faucet included with Quartz Top**

inside dimensions

10<sup>3</sup>/<sub>8</sub>" x 12<sup>3</sup>/<sub>8</sub>" x 6<sup>3</sup>/<sub>4</sub>"

26 cm x 31 cm x 17 cm

18 ga stainless steel



**DOOR & DRAWER SECURITY**

Add security to your casework with door and drawer locks, or lock and inside latch combination.

All locks on a cabinet are keyed alike unless otherwise requested.

**055** Door & Drawer Lock Option

**066** Door & Inside Latch Combo Option

**SIDE DESK**

Add more function to any size\* Fashion Finish casework with the addition of a wall/cabinet mounted desk.

**076FF** Desk with matching postform or quartz top

24" x 19" x 30"

61 cm x 48 cm x 76cm

\* Not available with quartz top with 8624





CLASSIC CLINTON

# CLASSIC CLINTON

♦ LAMINATE CASEWORK ♦

Classic Clinton Laminate Casework combines economical good looks with classic laminate construction. It is designed to provide physicians and clinics with cabinetry that meets storage needs and budgets.

Its easy clean laminate surface is scratch and scuff resistant and is available in popular woodgrain and solid finishes.

*Classic Clinton Laminate Casework features:*

- All laminated construction
- Easy-clean, durable, surfaces
- Soft-close door hinges
- 75 lbs. capacity drawers
- Choice of styles and sizes
- Easy installation
- Laminate countertops come with a 3 1/4" backsplash
- Optional stainless steel sink and faucet
- Optional postform countertops
- Choose from 5 laminates:

**Maple.**

Light and open, Maple blends soft hues of golden colors that brighten any decor.

**Dark Cherry.**

Use Dark Cherry's rich classic look to accent a variety of traditional styling.

**Gray**

Cool and contemporary, Gray is as versatile as it is trendy.

**Ashen Gray**

This gray has a warm feel and a rich hue that blends well with woods or other warm tones.

**Slate Gray**

A dark gray that goes well with bright accents, Slate Gray also works well as a contrasting countertop color.





Wall Cabinet with 2 Doors

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Easily mounts with hook and wall rail mounting system
- Shipped assembled for easy installation

8224	width	depth	height
	24"	12"	24"
	61 cm	30 cm	61 cm

Base Cabinet with 2 Doors & 2 Drawers

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Optional drawers may be added or repositioned to customize your layout
- Standard laminate, or optional postform laminate countertops
- Shown with optional sink and faucet
- Shipped assembled for easy installation

8024	width	depth	height
	24"	18"	35"
	61 cm	46 cm	89 cm



Wall Cabinet with 2 Doors

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Easily mounts with hook and wall rail mounting system
- Shipped assembled for easy installation

8236	width	depth	height
	36"	12"	24"
	91 cm	30 cm	61 cm

Base Cabinet with 2 Doors & 2 Drawers

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Optional drawers may be added or repositioned to customize your layout
- Standard laminate or optional postform laminate countertops
- Shown with optional sink and faucet
- Shipped assembled for easy installation

8036	width	depth	height
	36"	18"	35"
	91 cm	46 cm	89 cm



Wall Cabinet with 2 Doors

- 1 adjustable shelf
- Shipped assembled for easy installation
- Concealed, self-closing, adjustable, door hinges
- Easily mounts with hook rail mounting system

8242	width	depth	height
	42"	12"	24"
	107 cm	30 cm	61 cm

Base Cabinet with 3 Doors & 2 Drawers

- 2 adjustable shelves
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Optional drawers may be added or repositioned to customize your layout
- Standard laminate, or optional postform laminate countertops
- Shown with optional sink and faucet
- Shipped assembled for easy installation

8042	width	depth	height
	42"	18"	35"
	107 cm	46 cm	89 cm

Wall Cabinet with 2 Doors

- 1 adjustable shelf
- Shipped assembled for easy installation
- Concealed, adjustable, soft-close, door hinges
- Easily mounts with hook and rail mounting system

8248	width	depth	height
	48"	12"	24"
	122 cm	30 cm	61 cm

Base Cabinet with 4 Doors & 2 Drawers

- 2 adjustable shelves
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Optional drawers may be added or repositioned to customize your layout
- Standard laminate, or optional postform laminate countertops
- Shown with optional sink and faucet
- Shipped assembled for easy installation

8048	width	depth	height
	48"	18"	35"
	122 cm	46 cm	89 cm



Wall Cabinet with 3 Doors

- 2 adjustable shelves
- Shipped in 2 units for easy installation
- All holes are pre-drilled and all hardware is provided to connect the 2 cabinets
- Concealed, adjustable, soft-close, door hinges
- Easily mounts with hook and rail mounting system

8272	width	depth	height
	72"	12"	24"
	183 cm	30 cm	61 cm

Base Cabinet with 6 Doors & 3 Drawers

- 3 adjustable shelves
- Shipped assembled for easy installation
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Optional drawers may be added or repositioned to customize your layout
- Laminate or optional to postform countertops
- Shown with optional sink and faucet
- Shipped assembled for easy installation

8072	width	depth	height
	72"	18"	35"
	183 cm	46 cm	89 cm

Wall Cabinet with 1 Door


- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Wall cabinets shipped assembled for easy installation.
- Easy hook and wall rail mounting system

8224	width	depth	height
	24"	12"	24"
	61 cm	30 cm	61 cm

Base Cabinet with 2 Doors & Desk

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Standard laminate, or optional postform laminate countertops
- Shown with optional sink and faucet
- Base cabinet shipped assembled for easy installation. Desk easily attaches to cabinet

8024-99	width	depth	height
Overall	48"	18"	35"
	122 cm	46 cm	89 cm
Desk	24"	18"	30"
	61 cm	46 cm	76 cm



Wall Cabinet with 2 Doors

- 1 adjustable shelf
- Concealed, adjustable, soft-close, door hinges
- Wall cabinets shipped assembled for easy installation.
- Easy hook and wall rail mounting system

8242-9	width	depth	height
	42"	12"	24"
	107 cm	30 cm	61 cm

Base Cabinet with 3 Doors, 2 Drawers & Desk

- 2 adjustable shelves
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Optional drawers may be added or repositioned to customize your layout
- Standard laminate, or optional postform laminate countertops
- Shown with optional sink and faucet
- Shipped assembled for easy installation
- Base cabinet shipped assembled for easy installation. Desk easily attaches to cabinet

8042-99	width	depth	height
Overall	66"	18"	35"
	168 cm	46 cm	89 cm
Desk	24"	18"	30"
	61 cm	46 cm	76 cm

Wall Cabinet with 2 Doors

- 1 adjustable shelf
- Shipped assembled for easy installation
- Concealed, adjustable, soft-close, door hinges
- Easily mounts with hook and rail mounting system

8248	width	depth	height
	48"	12"	24"
	122 cm	30 cm	61 cm

Base Cabinet with 2 Doors, 5 Drawers & Desk

- 2 adjustable shelves
- Concealed, adjustable, soft-close, door hinges
- Adjustable, self-closing drawers hold up to 75 lbs.
- Standard laminate, or optional postform laminate countertops
- Shown with optional sink and faucet
- Base cabinets shipped assembled for easy installation. Desk easily attaches to cabinet

8048-99	width	depth	height
Overall	72"	18"	35"
	168 cm	46 cm	89 cm
Desk	24"	18"	30"
	61 cm	46 cm	76 cm





STANDARD LAMINATE COUNTERTOP

Proven in the field, Classic Clinton Laminate Countertops have been our standard for decades. Great looking and durable, laminate countertops are included with all Clinton Classic Casework unless a postform option is chosen. Classic Clinton Laminate Countertops can accent, or match your choice of laminate casework. A 3 1/4" backsplash is included with standard countertops


*NOTE: All base cabinets can accommodate an optional sink. Sinks on a 72" top may be positioned left, right or centered. 72" cabinet's tall doors MUST be positioned under sink to accommodate plumbing.*



OPTIONAL POSTFORM LAMINATE COUNTERTOP

Postform laminate countertops offer a distinctive one-piece countertop that eliminates the seam between countertop surface, backsplash and front edge. Postform countertops look great, are economical, durable and easy to clean.

*NOTE: Optional Postform countertops must be ordered separately, and are available only in the sizes indicated at right. All postform countertops without sinks are special order only. All 72" postform tops with sinks positioned other than centered are special order. 72" cabinet's tall doors MUST be positioned under sink to accommodate plumbing.*



MP	Maple
DC	Dark Cherry
GR	Gray
AS	Ashen Gray
SG	Slate Gray

24P	24" Top with sink & faucet
36P	36" Top with sink & faucet
42P	42" Top with sink & faucet
48P	48" Top with sink & faucet
72P	72" Top with sink & faucet



WC	White Cararra
BA	Black Alicante
ME	Meteorite

TOP MOUNT STAINLESS STEEL SINK AND FAUCET (LAMINATE OR POSTFORM)

Add more function to your cabinetry with a stylish, easy clean, top mounted, stainless steel sink and chrome faucet. Faucet is 14.5" high chrome with a 3" radius and dual, wing levers.

022 Sink & Faucet Option with Classic Laminate Top

Sink & Faucet included with Postform Top

inside sink dimensions
10" x 12" x 6 3/4"
25 cm x 30 cm x 17 cm
18 ga stainless steel



DOOR & DRAWER SECURITY

Add security to your casework with door and drawer locks, or lock and inside latch combination. All locks on a cabinet are keyed alike unless otherwise requested.

055 Door & Drawer Lock Option

066 Door & Inside Latch Combo Option

SIDE DESK

Add more function to any size Classic Clinton casework with the addition of a wall/cabinet mounted desk.

076 Desk

24" x 19" x 30"
61 cm x 48 cm x 76cm





Clinton Bedside Cabinets Treatment Cabinets and Mobile Carts are designed with the needs of the medical profession in mind. Built tough, with a wide range of sizes, models and a variety of storage configurations, Clinton bedside, treatment cabinets and carts are the perfect complement to any medical facility.

All models available in Clinton Classic laminates.  
Select models available in Fashion Finish.

Clinton Bedside and Treatment Cabinets feature:

- Easy-clean, all laminate construction, or laminate interior with fashion finish exteriors on select models
- Tough, chip-resistant edges on laminate cabinets and no front facing seams on Fashion Finish models
- Cabinets are assembled with steel fasteners, not glued wooden dowels
- 3-way adjustable, 110°, soft-close, door hinges
- Drawers with metal sides that hold up to 75 lbs.
- Recessed toe kick or durable caster bases
- Adjustable drawer fronts
- Mix color choices for base and top on select models
- Select models are available with a one-piece, gray, ABS molded top rimmed on all sides
- All laminated panels are California CARB II compliant
- 200 lbs. load capacity under normal use



8711



8711-A



8720



8720-A

Bedside Cabinet

- 1 self-closing drawer and door
- 1 adjustable shelf
- 2" high rimmed top on 3 sides (5 cm)
- Laminate drawer glides on nylon rollers
- Easy-clean, laminate surfaces
- Concealed, adjustable, soft-close, door hinges

8711	width	depth	height
	18 "	15 3/4"	28 3/4"
	44 cm	40 cm	73 cm

Available in Clinton Classic Laminate only.

Molded Top, Bedside Cabinet

- 1 self-closing drawer and door
- 1 adjustable shelf
- One-piece, vacuum-formed, gray, haircell textured top rimmed on all sides
- Laminate drawer glides on nylon rollers
- Easy-clean, laminate base
- Concealed, adjustable, soft-close, door hinges

8711-A	width	depth	height
	18"	16"	28 3/4"
	44 cm	41 cm	73 cm

Available in Clinton Classic Laminate only.

Mobile, Bedside Cabinet

- 1 self-closing drawer and door
- 1 adjustable shelf
- 2" high rimmed top on 3 sides (5 cm)
- Laminate drawer glides on nylon rollers
- Easy-clean, laminate surfaces
- Concealed, adjustable, soft-close, door hinges
- Mounted on 2" dual wheel, swivel casters

8720	width	depth	height
	18"	15 3/4"	29 1/4"
	43.81 cm	40 cm	74.3 cm

Available in Clinton Classic Laminate only.

Mobile, Molded Top, Bedside Cabinet

- 1 self-closing drawer and door
- 1 adjustable shelf
- One-piece, vacuum-formed, gray, haircell textured top rimmed on all sides
- Laminate drawer glides on nylon rollers
- Easy-clean, laminate base
- Concealed, adjustable, soft-close, door hinges
- Mounted on 2" dual wheel, swivel casters

8720-A	width	depth	height
	18"	16"	29 1/4"
	46 cm	41 cm	74 cm

Available in Clinton Classic Laminate only.



TREATMENT CABINETS

Mobile, Treatment Cabinet with 2 Doors and 1 Drawer

- 1 adjustable shelf
- Self-closing, smooth-glide, drawer with metal sides holds up to 75 lbs. (34 kg)
- All laminate base
- Concealed, adjustable, soft-close, door hinges
- All laminate top with 3 1/4" backsplash (8.25 cm)
- Mounted on 2" dual wheel, swivel casters

8921	width	depth	height
	24"	18 1/4"	35 1/4"
	61 cm	46 cm	90 cm

Available in Clinton Classic Laminate only.



8921

Mobile, Treatment Cabinet with 2 Doors and 2 Drawers

- 1 adjustable shelf
- Self-closing, smooth-glide, drawers with metal sides hold up to 75 lbs. (34 kg)
- All laminate base
- Concealed, adjustable, soft-close, door hinges
- All laminate top with 3 1/4" backsplash (8 cm)
- Mounted on 2" dual wheel, swivel casters

8922	width	depth	height
	24"	18 1/4"	35 1/4"
	61 cm	46 cm	90 cm

Available in Clinton Classic Laminate only.



8922

Mobile, Treatment Cabinet with 5 Drawers

- Self-closing, smooth-glide, drawers with metal sides hold up to 75 lbs. (34 kg)
- All laminate base
- All laminate top with 3 1/4" backsplash (8 cm)
- Mounted on 2" dual wheel, swivel casters

8950	width	depth	height
	24"	18 1/4"	35 1/4"
	61 cm	46 cm	90 cm

Available in Clinton Classic Laminate only.



8950



8805

Treatment Cabinet with 5 Drawers

- Smooth glide tracks for easy open and close
- Self-closing, smooth-glide, drawers with metal sides hold up to 75 lbs. (34 kg)
- All laminate base
- All laminate top with 3 1/4" backsplash (8 cm)

8805	width	depth	height
	24"	18 1/4"	35"
	61 cm	46 cm	89 cm

Available in Clinton Classic Laminate only.

Treatment Cabinet with 2 Doors & 1 Drawer

- 1 adjustable shelf
- Self-closing, smooth-glide, drawer with metal sides holds up to 75 lbs. (34 kg)
- Concealed, adjustable, soft-close, door hinges
- All laminate base
- All laminate top with 3 1/4" backsplash (8 cm)

8821	width	depth	height
	24"	18 1/4"	35"
	61 cm	46 cm	89 cm

Available in Clinton Classic Laminate only.



8821

Treatment Cabinet with 2 Doors & 2 Drawers

- 1 adjustable shelf
- Self-closing, smooth-glide, drawers with metal sides hold up to 75 lbs. (34 kg)
- Concealed, adjustable, soft-close, door hinges
- All laminate base
- All laminate top with 3 1/4" backsplash (8 cm)

8822	width	depth	height
	24"	18 1/4"	35"
	61 cm	46 cm	89 cm

Available in Clinton Classic Laminate only.



8822



TREATMENT CABINETS

Molded Top, Treatment Cabinet with 2 Doors & 1 Drawer

- Laminate interior construction
- Laminate or optional Fashion Finish exterior
- One-piece, vacuum-formed, gray haircell textured top rimmed on all sides,
- 1 adjustable shelf
- Self-closing, smooth-glide, drawer with metal sides holds up to 75 lbs. (34 kg)
- Concealed, adjustable, soft-close, door hinges

8821-A	width	depth	height
8821-AF	24 <sup>1</sup> / <sub>2</sub> " 62 cm	18 <sup>3</sup> / <sub>8</sub> " 47 cm	35" 89 cm



8821-A/Laminate



8821-AF/Fashion Finish



8921-A/Laminate



8921-AF/Fashion Finish

Moulded Top, Mobile, Treatment Cabinet with 2 Doors & 1 Drawer

- Laminate interior construction
- Laminate or optional Fashion Finish exterior
- One-piece, vacuum-formed, gray, haircell textured top rimmed on all sides
- 1 adjustable shelf
- Self-closing, smooth-glide, drawer with metal sides holds up to 75 lbs. (34 kg)
- Concealed, adjustable, soft-close, door hinges
- Clinton Classic Laminate cabinet is mounted on 2" dual wheel, swivel casters
- Fashion Finish cabinet is mounted on 3" dual wheel, locking, swivel casters

8921-A	width	depth	height
8921-AF	24 <sup>1</sup> / <sub>2</sub> " 62 cm	18 <sup>3</sup> / <sub>8</sub> " 47 cm	35" 89 cm

Molded Top, Treatment Cabinet with 2 Doors & 2 Drawers

- Laminate interior construction
- Laminate or optional Fashion Finish exterior
- One-piece, vacuum-formed, gray, haircell textured top rimmed on all sides
- 1 adjustable shelf
- Self-closing, smooth-glide, drawers with metal sides hold up to 75 lbs. (34 kg)
- Concealed, adjustable, soft-close, door hinges

8822-A	width	depth	height
8822-AF	24 <sup>1</sup> / <sub>2</sub> " 62 cm	18 <sup>3</sup> / <sub>8</sub> " 47 cm	35" 89 cm



8822-A/Laminate



8822-AF/Fashion Finish



8922-A/Laminate



8922-AF/Fashion Finish

Moulded Top, Mobile, Treatment Cabinet with 2 Doors & 2 Drawers

- Laminate interior construction
- Laminate or optional Fashion Finish exterior
- One-piece, vacuum formed, gray, haircell textured top rimmed on all sides
- 1 adjustable shelf
- Self-closing, smooth-glide, drawers with metal sides hold up to 75 lbs. (34 kg)
- Concealed, adjustable, soft-close, door hinges
- Clinton Classic Laminate cabinet is mounted on 2" dual wheel, swivel casters
- Fashion Finish cabinet is mounted on 3" dual wheel, locking, swivel casters

8922-A	width	depth	height
8922-AF	24 <sup>1</sup> / <sub>2</sub> " 62 cm	18 <sup>3</sup> / <sub>8</sub> " 47 cm	35" 89 cm

Moulded Top, Treatment Cabinet with 5 Drawers

- Laminate interior construction
- Laminate or optional Fashion Finish exterior
- One-piece, vacuum-formed, gray, haircell textured top rimmed on all sides
- Smooth glide tracks for easy open and close
- Self-closing, smooth-glide, drawers with metal sides hold up to 75 lbs. (34 kg)

8805-A	width	depth	height
8805-AF	24 <sup>1</sup> / <sub>2</sub> " 62 cm	18 <sup>3</sup> / <sub>8</sub> " 47 cm	35" 89 cm



8805-A/Laminate



8805-AF/Fashion Finish



8950-A/Laminate



8950-AF/Fashion Finish

Moulded Top, Mobile, Treatment Cabinet with 5 Drawers

- Laminate interior construction
- Laminate or optional Fashion Finish exterior
- One-piece, vacuum-formed, gray, haircell textured top rimmed on all sides
- Self-closing, smooth-glide, drawers with metal sides hold up to 75 lbs. (34 kg)
- Self-closing drawers with bumpers
- Clinton Classic Laminate cabinet is mounted on 2" dual wheel, swivel casters
- Fashion Finish cabinet is mounted on 3" dual wheel, locking, swivel casters

8950-A	width	depth	height
8950-AF	24 <sup>1</sup> / <sub>2</sub> " 62 cm	18 <sup>3</sup> / <sub>8</sub> " 47 cm	35" 89 cm

MOBILE CARTS

- Moulded Top, Mobile, Equipment Cart**
- Laminate, cabinet section with one, 4" drawer
  - One-piece, vacuum-formed, gray, haircell textured top rimmed on all sides
  - Gray, powder-coated, steel legs with four, dual wheel locking casters
  - Ships fully assembled

8900-A	width	depth	height
	24 <sup>1</sup> / <sub>2</sub> "	18 <sup>3</sup> / <sub>8</sub> "	30"
	62 cm	47 cm	76 cm



8900-A

- Laminate Top, Mobile, Equipment Cabinet**
- All laminate cabinet mounted on a stable, 5-leg black nylon base with dual wheel casters
  - 2 self-closing, drawers with bumpers
  - All laminate cabinet and top rimmed on 3 sides
  - Ships knocked-down

8920	width	depth	height
	22"	17"	33"
	56 cm	43 cm	84 cm



8920

- Moulded Top, Mobile, Equipment Cabinet**
- Laminate cabinet section mounted on a stable, 5-leg black nylon base with dual wheel casters
  - One-piece, vacuum-formed, gray, haircell textured top rimmed on all sides
  - Mounted on a stable 5-leg black nylon base
  - 2 self-closing, drawers with bumpers
  - Ships knocked-down

8923-A	width	depth	height
	24 <sup>1</sup> / <sub>2</sub> "	18 <sup>3</sup> / <sub>8</sub> "	33"
	62 cm	47 cm	84 cm



8923-A

- Optional Basket**
- Wire basket with white epoxy finish

089	width	depth	height
	17 <sup>1</sup> / <sub>2</sub> "	6 <sup>3</sup> / <sub>4</sub> "	5"
	44 cm	17 cm	13 cm







Put medical information and entry at your fingertips with Clinton's Tec-Cart™ Mobile Workstations. Tec-Carts™ feature large, easy-clean, moulded plastic work surfaces, pneumatic height adjustment for ease of use by multiple persons and large premium casters.

Clinton Tec-Carts make it easy for medical personnel to take versatile work stations anywhere their busy day takes them.

All Tec-Cart™ Mobile Medical Workstations feature:

- Dent and ding resistant, moulded plastic work surfaces
- Pneumatic, height adjustment
- Premium, 3" casters for easy mobility
- Cast aluminum, 5-leg base



**Compact, Tec-Cart™ Mobile Work Station with Drawer**

- Smooth rolling, premium, 3" locking casters
- Vacuum-formed, gray, textured, ABS plastic, work surface rimmed on all sides
- 5-leg, 24" cast aluminum base with bumper guards
- Easy-clean surfaces
- Pneumatic height adjustment for sitting or standing
- Gray laminate cabinet
- Easy-clean, all plastic drawer with gray laminate front and wire pull

9800	width	depth	height range
	18"	15 <sup>3</sup> / <sub>4</sub> "	28 <sup>1</sup> / <sub>2</sub> " – 39"
	46 cm	40 cm	72 cm – 99 cm

drawer	length	width	depth
	14 <sup>3</sup> / <sub>4</sub> "	10 <sup>1</sup> / <sub>2</sub> "	3"
	37 cm	27 cm	8 cm

**Large, Tec-Cart™ Mobile Work Station with Drawer**

- Smooth rolling, premium, 3" locking casters
- Vacuum-formed, gray, textured, ABS plastic, work surface rimmed on all sides
- 5-leg, 24" cast aluminum base with bumper guards
- Easy-clean surfaces
- Pneumatic height adjustment for sitting or standing
- Gray laminate cabinet
- Easy-clean, all plastic drawer with gray laminate front and wire pull

9810	width	depth	height range
	20"	20 <sup>1</sup> / <sub>2</sub> "	28 <sup>1</sup> / <sub>2</sub> " – 39"
	51 cm	52 cm	72 cm – 99 cm

drawer	length	width	depth
	14 <sup>3</sup> / <sub>4</sub> "	10 <sup>1</sup> / <sub>2</sub> "	3"
	37 cm	27 cm	8 cm

**Contour Tec-Cart™ Mobile Work Station**

- Smooth rolling, premium, 3" locking casters
- Curved contour, vacuum-formed, gray, textured, ABS plastic, work surface rimmed on the perimeter
- 5-leg, 24" cast aluminum base with bumper guards
- Easy-clean surface
- Pneumatic height adjustment for sitting or standing

9820	width	depth	height range
	26"	18 <sup>1</sup> / <sub>2</sub> "	30 <sup>1</sup> / <sub>2</sub> " – 40 <sup>5</sup> / <sub>8</sub> "
	66 cm	47 cm	77 cm – 103 cm



*Ask about Clinton's entire line of products.*

- Medical Tables
- Power Tables
- Seating & Medical Accessories
- Phlebotomy Equipment
- Pediatric Equipment
- Physical Therapy Equipment

525 East Market St., York, PA 17403

800-441-9131 • 717-848-3519

Toll-Free Fax 866-522-1233

Fax 717-843-5871

[www.clinton-ind.com](http://www.clinton-ind.com)