

2.7 SWC-015-041, SWC-015G-041, SWC-015-041G, SWC-015G-041G

Top Cabinet Size	4 ft ³ 25.5"Wx17.5"Dx15.5"H
Top Cabinet Size	10.75 ft ³ 25.5"Wx17.5"Dx41.5"H
Overall Size	30"W x 24"D x 74"H
Electrical	120VAC, 12Amp, 60Hz 8' cord, 15A provided
Heat Loss	Max 2360 Btu/hr

Standard Features

- ✓ Digital Temperature Control
- ✓ °F or °C Display
- ✓ Accuracy +/- 1°F (°C)
- ✓ Uniform Heating
- ✓ One adjustable shelf – Top Cabinet
- ✓ Two adjustable shelves – Bottom Cabinet
- ✓ Reverseable Door Hinges
- ✓ Door Lock
- ✓ Low-E Tempered Glass (glass door)

Options

- Seismic Anchoring
- Casters (2 with brakes, total 4)
- Leveling Feet (total 4)
- Shelf SS (additional)
- Left Side Door Hinge
- Audible Alarm (high & low temperature)
- Trim SS on One Side
- Trim SS on Two Sides
- Power Plug Hospital Grade

