

PRECISION MEDICAL®

ASPIRATORS



*Power & performance
when you need it!*



Portable Aspirator (PM66)

Choice of Power Source

The EasyGoVac models give you the option of a Lithium battery for faster re-charge times and no memory, a lead acid battery or an AC unit for lightweight non-portable use.

High Flow Rates

Even with filters, all models have high flow rates of between 20 and 26 LPM, with vacuum levels of 51 – 533mmHg.

Canister Options

All units are available with a dedicated 1200cc re-useable canister or a standard 800cc canister.

Built-in Battery Charger

All battery models have a built-in charger, eliminating the need to carry a heavy external charger which can also be lost. Everything is inside the unit, to charge it, all you need to do is plug it in.

Easy to Use

All models have an easy to read “power source indicator”, “battery test indicator” (Li model only) “low battery warning light” and a “fully charged indicator” to make operation by the caregiver and/or patient easier. A convenient carry bag is included with all battery models.



EasyGoVac
PM66Li-R

Part Number	Description
PM66AC	EasyGoVac Aspirator (AC only)
PM66Li	EasyGoVac Aspirator with Lithium Battery; includes Carry Bag
PM66Ped	Pediatric EasyGoVac Aspirator with Lithium Battery; includes Carry Bag
PM66S	EasyGoVac Aspirator with Sealed Lead Acid Battery; includes Carry Bag

When ordering a PM66 Aspirator add:
 -D for Disposable Canister
 -R for Reuseable Canister

Portable Aspirator (PM65)

It Performs ...

No need to worry about power when you need it. This high performance unit delivers a vacuum range of 2-21 inHg. The standard easy-to-read gauge allows the care giver to adjust the vacuum.

It Lasts...

The high performance battery allows 1.5 hours of run time, twice that of comparable units. A “Low Battery Light”, standard in this unit, informs the caregiver that the unit needs to be re-charged thus preventing interruptions in use.

It's Convenient ...

No need to buy a special canister. The Easy Go Vac uses any standard 800cc canister and is visible in the Carry Bag while unit is in operation. It is re-chargeable in its convenient, specially designed Carry Bag, which also stores the charger, Suction Tubing and optional Auto DC Power Cord.

It's Durable ...

The ABS plastic cabinet is durable, attractive, easy-to-clean and infection control friendly. There are no open areas on the unit to harbor bacteria and dirt.

It's Protected ...

A dual pump protection system features two safeguards. A float shut-off in the canister to protect overflow and an Inline Hydrophobic Filter to prevent back flow of fluid as well as airborne contaminants from entering the pump.

It's Dependable ...

The Easy Go Vac carries a full 2 year warranty on the unit and a 6 month warranty on the battery.



EasyGoVac
PM65

Part Number	Description
PM65	EasyGoVac Aspirator - Home Grade Charger
PM65HG	EasyGoVac Aspirator - Hospital Grade Charger

Stationary Aspirators (PM60, PM61 and PM63)

They are Powerful...

Our aspirators have vacuum ranges which are among the highest in the industry. The EasyVac (PM60), PowerVac (PM61) & PowerVac+ (PM63) aspirators have 1/5 horsepower performance.

They are Protected...

A dual protection system consists of a float shut-off in the canister lid and an Inline Hydrophobic Filter to prevent back flow of fluid as well as airborne contaminants from entering the pumps.

They are Versatile...

The EasyVac (PM60) has a protective canister bracket which not only holds the canister but also protects the gauge and adjustment knob from accidents. Safety is assured with the automatic resetting thermal limit switch which protects the pump from overheating.

The PM61 & PM63 have ABS plastic cabinets which are durable, attractive, easy to clean and infection control friendly. The PM63 has both Continuous and Intermittent Modes.

They are Safe...

Safety is assured with the automatic resetting thermal limit switch which protects the unit from overheating.



**EasyVac
PM60**



**PowerVac Aspirator
PM61**



**PowerVac+ Aspirator
PM63**

Part Number	Description
PM60	EasyVac Aspirator with Hospital Grade Power Cord (120VAC)
PM61	PowerVac Aspirator with Hospital Grade Power Cord & 2000 cc canister
PM63	PowerVac+ Aspirator with Hospital Grade Power Cord & 2000 cc canister

Aspirator Accessories

Canisters



- 502519-10 Canister (with Lid), 800 cc (Case of 10)
- 1906-12 Canister (with Lid), 1200 cc (Case of 12)
- 502687-12 Canister (with Lid), 2000 cc (Case of 12)
- 507276 Reusable Canister (autocalvable - 1500cc)
- 507522 1200cc Canister, reusable (PM66)

Tubing and Filters



- 1955 Patient Suction Tubing (6 ft)
- 1690 Inline Hydrophobic Filter
- 506047 Hydrophobic Filter (PM65)
- 502690 Inline Hydrophobic Filter & Tube Assembly
- 507754 Inline Filter (Bacteria) (PM66)
- 506851 Inline Filter (Hydrophobic) (PM66)

Power Cords



- 502954 Auto DC Power Cord (PM65)
- 507213 Auto DC Power Cord (PM66)
- 506961 AC Power Cord (PM66)

Roll Stands



- PM60-1 Roll Stand
(for use with any PM aspirators or compressors)
- PM60-2 Roll Stand with canister bracket
(for use with any PM aspirators)

Inlet Adaptor Kit



- 507521 800cc Inlet Adapter Kit (PM66 only)

Carry Bags



- 507263 Carry Bag (PM66)
- 502956 Carry Bag (PM65)

Battery Charger/Power Supply



- 503974 Battery Charger/Power Supply (PM65 only)

Canister Bracket



- 502540 Canister Bracket (PM61 & PM63 only)